

TEST CHANGE

Vaginitis Panel by TMA

3002581, VPAN TMA

Specimen Requirements:

Patient Preparation:	Patient must be 14 years of age or older.
Collect:	Vaginal specimen collected with pink swab from Aptima MultiTest Swab Collection kit (ARUP supply #65761 single collection kit #55224 PK/50 or #55224 pack of 50 collection kits 55229 PK/10) available online through eSupply using ARUP Connect(TM) or contact Client Services at (800-)522-2787.
Specimen Preparation:	Place swab in MultiTest Swab Specimen Transport Tube, break shaft at scoreline, then recap tube.
Transport Temperature:	Refrigerated.
Unacceptable Conditions:	Specimens in any transport media other than indicated above. Specimen in MultiTest swab transport media without a swab.

Remarks:

Stability: Ambient: 30 days; Refrigerated: 30 days; Frozen: 90 days

Methodology: Qualitative Transcription-Mediated Amplification [\(TMA\)](#)

Note:

CPT Codes: 81513; 87481 x2; 87661

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

See report

Reference Interval:

Test Number	Components	Reference Interval
	Bacterial Vaginosis by TMA	Negative
		Negative
	Candida glabrata by TMA	Negative
		Negative
	Candida species (other) by TMA	Negative
		Negative
	Trichomonas vaginalis by TMA	Negative
		Negative

