
3002850

HLA Antibody Screen, Class I and Class II

HLAABSCN

Specimen Required: Collect: Plain red.

Specimen Preparation: Transfer 5 mL serum to ARUP Standard Transport Tubes. (Min. 2 mL)

New York State Clients: Transfer 7 mL serum to ARUP Standard Transport Tubes. (Min. 3 mL)

Storage/Transport Temperature: Refrigerated

Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 month; Frozen: 2 years