

HOTLINE: Effective August 15, 2022

**3002887 Autoimmune Neurologic Disease Reflexive Panel, CSF**
**NEURORCSF**

**Methodology:** Semi-Quantitative **Cell-Based** Indirect Fluorescent Antibody/Qualitative Immunoblot/Quantitative Radioimmunoassay/Semi-quantitative Enzyme-Linked Immunosorbent Assay

**Reference Interval:**

Test Number	Components	Reference Interval
2005164	N-methyl-D-Aspartate Receptor Antibody, IgG <b>CBA-IFA</b> , CSF with Reflex to Titer	Effective May 21, 2012 < 1:1
3001257	Alpha-Amino-3-hydroxy-5-methyl-4-isoxazolepropionic Acid (AMPA) Receptor Antibody, IgG by <b>CBA-IFA</b> with Reflex to Titer, CSF	Less than 1:1
3001267	Gamma Aminobutyric Acid Receptor, Type B (GABA-BR) Antibody, IgG by <b>CBA-IFA</b> with Reflex to Titer, CSF	Less than 1:1
3001986	Contactin-Associated Protein-2 Antibody, IgG <b>CBA-IFA</b> with Reflex to Titer, CSF	Less than 1:1
3001387	Voltage-Gated Potassium Channel (VGKC) Antibody, CSF	
		Negative 0.0-1.1 pmol/L
		Positive 1.2 pmol/L or greater
2010841	Paraneoplastic Antibodies (PCCA/ANNA) by IFA with Reflex to Titer and Immunoblot, CSF	
3001992	Leucine-Rich, Glioma-Inactivated Protein 1 Antibody, IgG <b>CBA-IFA</b> with Reflex to Titer, CSF	Less than 1:1
3002257	CV2.1 Screen by <b>CBA-IFA</b> with Reflex to Titer, CSF	Less than 1:1
3002788	Glutamic Acid Decarboxylase Antibody, CSF	0.0-5.0 IU/mL
3002886	SOX1 Antibody, IgG by Immunoblot, CSF	Negative
3004512	Dipeptidyl Aminopeptidase-Like Protein 6 (DPPX) Antibody, IgG by <b>CBA-IFA</b> With Reflex to Titer, CSF	Less than 1:1
3004510	Amphiphysin Antibody IgG, CSF	Negative