

HOTLINE: Effective August 15, 2022

3002929 Paraneoplastic Reflexive Panel

PNS PAN2

Methodology: Semi-Quantitative Cell-Based Indirect Fluorescent Antibody/Qualitative Immunoblot

Reference Interval:

Test Number	Components	Reference Interval		
2013956	CV2.1 Screen by CBA-IFA with Reflex to Titer	Less than 1:10		
2007961	PCCA/ANNA by IFA with Reflex to Titer and Immunoblot	Effective August 17, 2020		
		Test Number	Components	Reference Interval
			Purkinje Cell/Neuronal Nuclear IgG Scrn	None Detected
			Neuronal Nuclear Antibody (ANNA) IFA Titer, IgG	Less than 1:10
			Purkinje Cell Antibody, Titer	Less than 1:10
		3002917	Neuronal Nuclear Antibodies (Hu, Ri, Yo, Tr/DNER) IgG by Immunoblot, Serum	Refer to report
2008893	Amphiphysin Antibody, IgG	Negative		
3002885	SOX1 Antibody, IgG by Immunoblot, Serum	Negative		