

HOTLINE: Effective February 1, 2021

3003039

Cyanide, Whole Blood

CYANI WB

Specimen Required: Collect: Gray top tube (Sodium Fluoride / Potassium Oxalate)

Specimen Preparation: 1 mL whole blood. (Min: 0.4 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.

Stability (collection to initiation of testing): Ambient: Undetermined; Refrigerated: 1 week; Frozen: 3 months