

TEST CHANGE

Kratom (Mitragynine) - Screen ~~W~~with Reflex to Confirmation/~~Quantitation~~, Urine
3003311, KRATOM U

Specimen Requirements:

Patient Preparation:

Collect: Urine.

Specimen Preparation: Transfer 10 mL urine to ARUP ~~standard transport tubes~~~~Standard Transport Tubes~~. (Min: 1 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen. Also acceptable: Refrigerated.

Unacceptable Conditions: Urine from preservative tube.

Remarks:

Stability: Ambient: 72 hours; Refrigerated: 2 weeks; Frozen: 6 months

Methodology: Qualitative Immunoassay/~~Qualitative Liquid Chromatography-Tandem Mass Spectrometry~~

Performed: Varies

Reported: 6-9 days

Note: If screen is positive, then confirmation will be added. Additional charges apply.

CPT Codes: 80307 if reflexed, add 80323

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

By report