

3003634

**Capillary Malformation-Arteriovenous Malformation (CM-AVM) Panel,
Sequencing and Deletion/Duplication**

CMAVM NGS

Methodology: Massively Parallel Sequencing
Performed: Varies
Reported: 3 weeks

Specimen Required: Collect: Lavender (EDTA) or yellow (ACD solution A or B).
Specimen Preparation: Transport 3 mL whole blood. (Min: 3 mL)
Storage/Transport Temperature: Refrigerated
Unacceptable Conditions: Serum or plasma; grossly hemolyzed or frozen specimens; saliva, buccal brush, or swab; FFPE tissue.
Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: Unacceptable

Interpretive Data:
Refer to report

Note: *EPHB4* (NM_004444); *RASA1* (NM_002890)