

HOTLINE: Effective May 17, 2021

---

**New Test**      **3003676**      **Vedolizumab Quantitation with Antibodies, Serum**      **VEDOL AB**

**Methodology:**      Quantitative Liquid Chromatography/Mass Spectrometry (LC-MS/MS) /Electrochemiluminescent Immunoassay  
**Performed:**      Varies  
**Reported:**      7-17 days

**Specimen Required:** Patient Prep: 12 hours prior to specimen collection discontinue multivitamins or dietary supplements containing biotin (vitamin B7), commonly found in hair, skin, and nail supplements. Nivolumab (Opdivo) must be discontinued at least 4 weeks prior to testing.  
Collect: Plain red. Also acceptable: Serum separator tube (SST). Collect immediately before next scheduled dose (trough specimen).  
Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1.5 mL serum to an ARUP Standard Transport Tube. (Min: 0.75 mL)  
**Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.**  
Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen.  
Stability (collection to initiation of testing): Ambient: Unacceptable; Refrigerated: 28 days; Frozen: 28 days

**Reference Interval:**      By Report

**CPT Code(s):**      80280; 82397

New York DOH Approved.

**HOTLINE NOTE:** Refer to the Test Mix Addendum for interface build information.