

HOTLINE: Effective May 17, 2021

---

**New Test**      **3003726**      **Propoxyphene and Metabolite, Urine**      **PROPOX U**

**Methodology:** Quantitative Gas Chromatography-Mass Spectrometry (GC-MS)  
**Performed:** Varies  
**Reported:** 8 – 11 days

**Specimen Required:** Collect: Urine  
Specimen Preparation: Transfer 2 mL urine to an ARUP Standard Transport Tube. (Min: 0.7 mL)  
**Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.**  
Storage/Transport Temperature: Refrigerated. Also acceptable: Frozen  
Stability (collection to initiation of testing): Ambient: 72 hours; Refrigerated: 1 week; Frozen: 3 months

**Reference Interval:** By report

**Note:** Amitriptyline is a known interference.

**CPT Code(s):** 80367 (Alt code: G0480)

New York DOH Approved.

**HOTLINE NOTE:** Refer to the Test Mix Addendum for interface build information.