

TEST CHANGE

JAK2 (V617F) Mutation by ddPCR, Quantitative

3003751, JAK2V617FQ

Specimen Requirements:

Patient Preparation:

Collect: Whole blood or bone marrow in lavender (EDTA) **or pink (K2EDTA)**. Also acceptable: Whole blood in green (sodium heparin).

Specimen Preparation: Whole Blood: Do not freeze. Transport 5 mL whole blood. (Min: 1 mL)
Bone Marrow: Do not freeze. Transport 3 mL bone marrow. (Min: 1 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma, serum, FFPE tissue blocks/slides, or fresh or frozen tissue. Specimens collected in anticoagulants other than EDTA or sodium heparin. Clotted or grossly hemolyzed specimens.

Remarks:

Stability: Refrigerated: 7 days; Frozen: Unacceptable

Methodology: Droplet Digital PCR (ddPCR)

Note:

CPT Codes: 81270

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

Reference Interval: