

HOTLINE: Effective May 16, 2022

3003801 JAK2 (V617F) Mutation by ddPCR, Qualitative with Reflex to JAK2 Exon 12 **Mutation Analysis by PCR**

PV RFX

Specimen Required: Collect: Whole blood or bone marrow: Lavender (EDTA).

Specimen Preparation: Whole Blood: Do not freeze. Transport 5 mL whole blood. (Min: 1 mL)

Bone Marrow: Do not freeze. Transport 3 mL bone marrow. (Min: 1 mL)

Storage/Transport Temperature: Refrigerated.

Unacceptable Conditions: Plasma, serum, FFPE tissue blocks/slides, or fresh or frozen tissue. Specimens collected in anticoagulants

other than EDTA. Clotted or grossly hemolyzed specimens.

Stability (collection to initiation of testing): Refrigerated: 4 days; Frozen: Unacceptable