

HOTLINE: Effective July 5, 2022

3003913

Orthopedic Metals Panel (Chromium, Cobalt, Titanium)

ORTHO PAN

Specimen Required: Collect: Body fluid.

Specimen Preparation: Transfer 5 mL body fluid to a trace element-free transport tube (ARUP supply #43116) or acid-washed transfer vial (ARUP supply #54350) available online through eSupply using ARUP Connect™ or contact ARUP Client Services at (800) 522-2787. (Min: 2.2 mL)

Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Storage/Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Stability (collection to initiation of testing): Ambient: Undetermined; Refrigerated: Undetermined; Frozen: Undetermined