

HOTLINE: Effective September 14, 2021

---

<b>New Test</b> Available Now	<b>3004244</b>	<b>BAP1 by Immunohistochemistry</b>	<b>BAP1 IHC</b>
----------------------------------	----------------	-------------------------------------	-----------------

---

**Methodology:** Immunohistochemistry  
**Performed:** Mon-Fri  
**Reported:** 1-3 days

**Specimen Required:** Collect: Tissue.

Specimen Preparation: Formalin fix (10 percent neutral buffered formalin) and paraffin embed specimen (cells must be prepared into a cellblock). Protect paraffin block and/or slides from excessive heat. Transport tissue block or 5 unstained (3- to 5-micron thick sections), positively charged slides in a tissue transport kit (recommended but not required), (ARUP supply #47808) available online through eSupply using ARUP Connect or contact ARUP Client Services at (800) 522-2787. (Min: 2 slides) If sending precut slides, do not oven bake.

Storage/Transport Temperature: Room temperature. Also acceptable: Refrigerated. Ship in cooled container during summer months.

Unacceptable Conditions: Specimens submitted with non-representative tissue type. Depleted specimens.

Stability (collection to initiation of testing): Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

**Interpretive Data:**

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

**Note:** This test is performed as a stain and return (technical) service only.

**CPT Code(s):** 88342

New York DOH Approved.

**HOTLINE NOTE:** Refer to the Test Mix Addendum for interface build information.