

TEST CHANGE

Leukotriene E4, 24-Hour Urine 3004792, LTE 24 URN	
Specimen Requirements:	
Patient Preparation:	Patients taking 5-lipoxygenase inhibitor <u>z</u> Zileuton (/Zyflo) may have decreased concentrations of leukotriene E4 (LTE4). If possible, discontinue 48 hours prior to collection.
Collect:	24-hour urine. Refrigerate during collection.
Specimen Preparation:	From a well-mixed 24-hour collection transfer 5 mL urine to ARUP <u>standard transport tubes</u> <u>Standard Transport Tubes</u> (Min: 2 mL). Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.
Transport Temperature:	Frozen. Also acceptable: Refrigerated
Unacceptable Conditions:	
Remarks:	Collection duration and urine volume must be provided for testing. Total collection volume must be greater than 500 mL.
Stability:	Ambient: 24 hours; Refrigerated: 1 week; Frozen: <u>28 days</u> 1- month
Methodology:	Quantitative Colorimetry/Liquid Chromatography-Tandem Mass Spectrometry
Performed:	Varies
Reported:	3- <u>12</u> 9 days
Note:	
CPT Codes:	82542
New York DOH Approval Status:	This test is New York DOH approved.
Interpretive Data:	
Reference Interval:	
By Report	