

## NEW TEST

[Click for Pricing](#)

### Phospho-Tau/Total-Tau/A Beta42, CSF

3016444, PHOSPHO T

#### Specimen Requirements:

##### Patient Preparation:

Collect: CSF

Specimen Preparation: Transfer 2 mL CSF to an ARUP standard transport tube. (Min: 0.5 mL).

Transport Temperature: Frozen

Unacceptable Conditions: Specimens with cell count greater than 500 erythrocytes/mm<sup>3</sup>.

##### Remarks:

Stability: Ambient: 72 hours, Refrigerated: 21 days, Frozen: 4 months

Methodology: Enzyme-Linked Immunosorbent Assay (ELISA)

Performed: Varies

Reported: 7-19 days

##### Note:

CPT Codes: 83520 x3

New York DOH Approval Status: This test is New York DOH approved.

##### Interpretive Data:

##### Reference Interval:

Test Number	Components	Reference Interval
-------------	------------	--------------------

**HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.**