

NEW TEST

[Click for Pricing](#)

CSF, Glucose

3016614, CSF GLU

Specimen Requirements:

Patient Preparation:

Collect: CSF. Also acceptable: CSF collected in plain red.

Specimen Preparation: Centrifuge to remove cellular material. Transfer 1 mL CSF to an ARUP standard transport tube. (Min. 0.2 mL)

Transport Temperature: Refrigerated

Unacceptable Conditions: Specimen types other than those listed.

Remarks:

Stability:

Methodology: Quantitative Enzymatic Assay

Performed: Sun-Sat

Reported: Within 24 hours

Note:

CPT Codes: 82945

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval
	CSF, Glucose	50-80 mg/dL

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.