

NEW TEST

[Click for Pricing](#)

Anti-Dense Fine Speckled 70 (DFS70) Antibody by ELISA, Serum
3016769, DFS70 AB

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serum separator tube (SST).

Specimen Preparation: Transfer 1 mL serum to an ARUP standard transport tube. (Min: 0.5 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Frozen. Also acceptable: Refrigerated.

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: 2 weeks; Frozen: 6 months

Methodology:

Performed: Varies

Reported: 7-14 days

Note:

CPT Codes: 86038

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.