

NEW TEST

Click for Pricing

Pepsin A, Respiratory 3016813, PEPSIN		
Specimen Requirements:		
Patient Preparation:		
Collect:	Bronchial wash, broncheoalveolar lavage (BAL), or tracheal aspirate.	
Specimen Preparation:	Transfer 2 mL respiratory specimen to an ARUP standard transport tube and freeze immediately. (Min: 0.5 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.	
Transport Temperature:	Frozen	
Unacceptable Conditions:		
Remarks:		
Stability:	Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 3 months	
Methodology:		
Performed:	Varies	
Reported:	5-10 days	
Note:		
CPT Codes:	83986, 84157, 83516	
New York DOH Approval Status:	Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.	
Interpretive Data:		
Reference Interval:		
Test Components Number		Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.