

**NEW TEST – Available Now**

[Click for Pricing](#)

**CD30 for Bone Marrow Specimens by Immunohistochemistry**

3016879, CD30BM IHC

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Tissue or cells.

**Specimen Preparation:** Formalin fix (10 percent neutral buffered formalin) and paraffin embed specimen (cells must be prepared into a cellblock). Protect paraffin block and/or slides from excessive heat. Transport tissue block or 5 unstained (3- to 5-micron thick sections), positively charged slides in a tissue transport kit (ARUP supply #47808). Available online through eSupply using ARUP Connector contact ARUP Client Services at 800-522-2787. (Min: 2 slides). If sending precut slides, do not oven bake.

**Transport Temperature:**

**Unacceptable Conditions:** Tissue or cells not processed and placed in a paraffin block; serum, blood, or other body fluids; tissue not mounted on positively charged slides.

**Remarks:**

**Stability:** Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

**Methodology:** Qualitative Immunohistochemistry (IHC)

**Performed:** Mon-Fri

**Reported:** 1-3 days

**Note:** This test is performed as a stain and return (technical) service only.

**CPT Codes:** 88342

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

**Reference Interval:**

Test	Components	Reference Interval
------	------------	--------------------

Number	
--------	--

**HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.**