

NEW TEST

Click for Pricing

Vitamin B7, Serum or Plasma

3016932, VITA B7

Specimen Requirements:

Patient Preparation:

Collect: Plain red or lavender (K2EDTA)

Specimen Preparation: Transfer 1 mL serum or plasma to an ARUP standard transport

tube (Min: 0.5 mL). Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Effective Date: October 2, 2023

Transport Temperature: Frozen

Unacceptable Conditions: Separator tubes

Remarks:

Stability: Ambient: 72 hours; Refrigerated: 72 hours; Frozen: 2 weeks

Methodology: Quantitative Liquid Chromatography-Tandem Mass

Spectrometry

Performed: Varies

Reported: 3-10 days

Note:

CPT Codes: 84591

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Components Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.