

NEW TEST – Available Now

Click for Pricing

NUDT15 Genotyping

3017373, NUDT15GENO

Specimen Requirements:

Patient Preparation:

Collect: Lavender (EDTA), pink (K2EDTA), or yellow (ACD solution A or

B).

Specimen Preparation: Transport 3 mL whole blood. (Min: 1 mL)

Transport Temperature: Refrigerated

Unacceptable Conditions: Plasma or serum; specimens collected in sodium heparin or

lithium heparin; frozen specimens in glass collection tubes;

Effective Date: May 20, 2024

frozen yellow (ACD solution A or B)

Remarks:

Stability: Room temperature: 72 hours; Refrigerated: 1 week; Frozen: 1

month

Methodology: Polymerase Chain Reaction (PCR)/Fluorescence Monitoring

Performed: Varies

Reported: 5-10 days

Note: Whole blood is the preferred specimen. Saliva samples that

yield inadequate DNA quality and/or quantity will be reported as inconclusive if test performance does not meet laboratory-

determined criteria for reporting.

CPT Codes: 81306

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

Counseling and informed consent are recommended for genetic testing. Consent forms are

available online.

Reference Interval:

By report



Effective Date: May 20, 2024

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.