

**NEW TEST**

[Click for Pricing](#)

**Platelet Antibody Identification Panel**

3017737, PAB PAN

**Specimen Requirements:**

**Patient Preparation:**

**Collect:** Plain red. Also acceptable: Serum separator tube (SST).

**Specimen Preparation:** Transfer 5 mL serum to an ARUP standard transport tube. (Min: 3 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

**Transport Temperature:** Frozen. Also acceptable: Refrigerated.

**Unacceptable Conditions:**

**Remarks:**

**Stability:** Ambient: Unacceptable; Refrigerated: 1 week; Frozen: 2 weeks

**Methodology:** Platelet Antibody Bead Array (PABA)/Flow Cytometry

**Performed:** Varies

**Reported:** 7-10 days

**Note:** Panel includes testing for antibodies against HPA-1 a/b, HPA-2a/b, HPA-3a/b, HPA-4b, HPA-5a/b, GPIIb/IIIa, GPIa/IIa, GP1b/IX, GPIV and class I HLA.

**CPT Codes:** 86022

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

**Reference Interval:**

Test Number	Components	Reference Interval

**HOTLINE NOTE:** Refer to the Hotline Test Mix for interface build information.