

**NEW TEST**

[Click for Pricing](#)

**Tobramycin Level, Trough, Serum**

3018762, TOBRA TR

**Specimen Requirements:**

**Patient Preparation:** Collect specimen 30 minutes before next dose.

**Collect:** Serum separator tube (SST). Also acceptable: Plain red

**Specimen Preparation:** Separate from cells ASAP or within 2 hours of collection. Transfer 0.5 mL serum to an ARUP standard transport tube. (Min: 0.25 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

**Transport Temperature:** Frozen. Also acceptable: Refrigerated

**Unacceptable Conditions:**

**Remarks:**

**Stability:** Ambient: 72 hours; Refrigerated: 1 week; Frozen: 28 days

**Methodology:** Quantitative Enzyme Immunoassay (EIA)

**Performed:** Varies

**Reported:** 3-5 days

**Note:**

**CPT Codes:** 80200

**New York DOH Approval Status:** This test is New York DOH approved.

**Interpretive Data:**

**Reference Interval:**

Test Number	Components	Reference Interval

**HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.**