

NEW TEST

[Click for Pricing](#)

Dermatomyositis Autoantibody Panel

3018870, DERM PAN2

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST)

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer one 1 mL serum aliquots to ARUP standard transport tubes. (Min: 0.5 mL/aliquot)

Transport Temperature: Refrigerated

Unacceptable Conditions: Hemolyzed, hyperlipemic, icteric, heat-treated, or contaminated specimens

Remarks:

Stability: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year

Methodology: Qualitative Immunoprecipitation / Qualitative Immunoblot / Semi-Quantitative Indirect Fluorescent Antibody (IFA)

Performed: Sun-Sat

Reported: 7-18 days

Note: Antibodies: Mi-2, P155/140, SAE1, MDA5, NXP2, TIF1-gamma, ANA

CPT Codes: 83516 x2; 84182 x4; 86039

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.