

NEW TEST

[Click for Pricing](#)

Hydroxychloroquine, Blood

3018999, HCQ

Specimen Requirements:

Patient Preparation:

Collect: Lavander (K2EDTA)

Specimen Preparation: Transport 1 mL whole blood. (Min: 0.2 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 month; Refrigerated: 1 month; Frozen: 1 month

Methodology: Quantitative High Performance Liquid Chromatography-Tandem Mass Spectrometry

Performed: Varies

Reported: 6-10 days

Note:

CPT Codes: 80220

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.