

NEW TEST

Click for Pricing	
Lithium, Serum/Plasma	
3019008, LITHIUM SP	
Specimen Requirements:	
Patient Preparation:	
Collect:	Plain red, lavender (K2EDTA), or pink (K2EDTA).
Specimen Preparation:	Separate from cells ASAP or within 2 hours of collection. Transfer 2 mL serum or plasma to an acid-washed transfer vial (ARUP supply #54350) available online through eSupply using ARUP Connect (TM) or contact ARUP Client Services at (800) 522-2787. (Min: 0.4 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.
Transport Temperature:	Frozen. Also acceptable: Refrigerated or room temperature.
Unacceptable Conditions:	Light green top tube (lithium heparin).
Remarks:	
Stability:	Ambient: 2 weeks; Refrigerated: 2 weeks; Frozen: 1 year
Methodology:	Quantitative Inductively Coupled Plasma-Mass Spectrometry (ICP-MS)
Performed:	Varies
Reported:	7-11 days
Note:	
CPT Codes:	80178
New York DOH Approval Status:	This test is New York DOH approved.
Interpretive Data:	
Reference Interval:	
Test Components Number	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.