

NEW TEST

[Click for Pricing](#)

Digoxin, Serum or Plasma

3019009, DIGOXI SP

Specimen Requirements:

Patient Preparation: Collect specimen no sooner than 6 hours after dose.

Collect: Serum separator tube (SST) or gray (sodium fluoride/potassium oxalate).

Specimen Preparation: Separate from cells ASAP or within 2 hours of collection. Transfer 1 mL serum or plasma to an ARUP standard transport tube. (Min: 0.4 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 28 days; Refrigerated: 28 days; Frozen: 1 year

Methodology: Quantitative High Performance Liquid Chromatography-Tandem Mass Spectrometry

Performed: Varies

Reported: 12-16 days

Note:

CPT Codes: 80162

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.