

NEW TEST

Click for Pricing

ABPA With Quantitative A. Fumigatus IgG

3019329, ABPA PAN

Specimen Requirements:

Patient Preparation:

Collect: Plain red or serum separator tube.

Specimen Preparation: Separate serum from cells within 2 hours of collection.

Transfer 0.5 mL serum to an ARUP standard transport tube.

Effective Date: April 21, 2025

(Min: 0.3 mL)

Transport Temperature: Refrigerated.

Unacceptable Conditions: Postmortem samples and NY samples (non-New York

approved)

Remarks:

Stability: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 year

Methodology: Quantitative ImmunoCAP Fluorescent Enzyme Immunoassay

Performed: Sun-Sat

Reported: 1-8 days

Note:

CPT Codes: 82785; 86003; 86317

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Components Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.