

TEST CHANGE

von Willebrand Factor (VWF) GPIbM Activity

3019671, VWF_GPIBM

Specimen Requirements:

Patient Preparation:

Collect: Light blue (sodium citrate). [Refer to Specimen Special Collection and Handling Hemostasis/Thrombosis Specimens guide located at https://aruplab.com for hemostasis/thrombosis specimen handling guidelines./Specimen-Handling/SpecialSpecimenCollection/Hemostasis-Thrombosis.pdf](https://aruplab.com/hemostasis/thrombosis-specimen-handling-guidelines/specimen-handling/special-specimen-collection/hemostasis-thrombosis.pdf)

Specimen Preparation: Transfer 1.5 mL platelet-poor plasma to an ARUP standard transport tube. (Min: 1.0 mL)

Transport Temperature: CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.

Unacceptable Conditions: Serum, EDTA plasma, clotted, or hemolyzed specimens

Remarks:

Stability: Ambient: 4 hours; Refrigerated: Unacceptable; Frozen: 3 months

Methodology: Quantitative Immunospectrophotometry

Note:

CPT Codes: 85397

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval	
	von Willebrand Factor, Activity (GPIbM)		
		Age	Reference Intervals (%)
		0-6 years	51-215
		7-9 years	52-176
		10-11 years	60-195
		12-13 years	50-184
		14-15 years	50-203
		16-17 years	49-204

		18 years and older	51-215
--	--	--------------------	--------