

NEW TEST

Click for Pricing

Lactate, CSF

3019672, LACT CSF

Specimen Requirements:

Patient Preparation:

Collect: CSF

Specimen Preparation: Centrifuge to remove cellular material. Transfer 1 mL CSF to an

ARUP standard transport tube. (Min: 0.5 mL) Test is not performed at ARUP; separate specimens must be submitted

Effective Date: July 21, 2025

when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN

Unacceptable Conditions:

Remarks:

Stability: Ambient: 3 hours; Refrigerated: 24 hours; Frozen: 2 months

Methodology: Quantitative Colorimetry

Performed: Varies

Reported: 7-11 days

Note:

CPT Codes: 83605

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Components Reference Interval

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.