

NEW TEST – Available Now

[Click for Pricing](#)

GI Bacterial Panel by PCR With Reflex to GI Expanded Bacterial Panel by PCR

3020663, GI BACTRFX

Specimen Requirements:

Patient Preparation:

Collect: Stool in liquid Cary-Blair transport media

Specimen Preparation: Transfer stool to enteric transport media (Cary-Blair) (ARUP supply #29799) available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at 800-522-2787. Refrigerate after transfer.

Transport Temperature: Refrigerated

Unacceptable Conditions: Unpreserved stool, stool in media other than liquid Cary-Blair.

Remarks:

Stability: Ambient: 4 days; Refrigerated: 7 days; Frozen: Unacceptable

Methodology: Qualitative Polymerase Chain Reaction (PCR)

Note:

If GI Bacterial Panel by PCR results are negative for all components, then GI Expanded Bacterial Panel by PCR will be added. Additional charges apply.

CPT Codes: 87505; if reflexed, add 87505

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Refer to report

Test Number	Components	Reference Interval
	Campylobacter by PCR	Not Detected
	Salmonella by PCR	Not Detected
	Shiga toxin-producing E. coli (STEC) PCR	Not Detected
	Shigella/Enteroinvasive E. coli by PCR	Not Detected

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.