

NEW TEST

[Click for Pricing](#)

Matrix Metalloproteinase-9 (MMP-9), Serum

3020899, MMP9 S

Specimen Requirements:

Patient Preparation:

Collect: Serum separator tube (SST)

Specimen Preparation: Transfer 1 mL serum to an ARUP standard transport tube and freeze immediately. (Min: 0.21 mL)
Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: CRITICAL FROZEN

Unacceptable Conditions:

Remarks:

Stability: Ambient: Unacceptable; Refrigerated: Unacceptable; Frozen: 1 month

Methodology: Enzyme-Linked Immunosorbent Assay (ELISA)

Note:

CPT Codes: 83520

New York DOH Approval Status: Specimens from New York clients will be sent out to a New York DOH approved laboratory, if possible.

Interpretive Data:

Reference Interval:

Refer to report

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.