

NEW TEST

[Click for Pricing](#)

MYD88 L265P Mutation Detection, Tissue, by PCR, Quantitative

3021230, MYD88TISS

Specimen Requirements:

Patient Preparation:

Collect: Fresh tissue or formalin-fixed tumor tissue.

Specimen Preparation: FFPE Tumor Tissue : Formalin-fixed (10 percent neutral buffered formalin) and paraffin-embedded tissue. Protect from excessive heat. Transport tissue in a tissue transport kit (ARUP Supply #47808) available online through eSupply using ARUP Connect(TM) or contact ARUP Client Services at 800-522-2787.

Transport Temperature: FFPE Tumor Tissue : Room temperature. Also acceptable: Refrigerated. Ship in cooled container during summer months.

Unacceptable Conditions: FFPE Tumor Tissue: Specimens fixed in any fixative other than 10 percent neutral buffered formalin.

Remarks:

If multiple specimens (blocks or slides) are sent to ARUP, they must be accompanied by one of the following: an order comment indicating that the ARUP pathologist should choose the specimen most appropriate for testing (e.g., "Choose best block"), or individual orders for each sample submitted. A Pathologist Block Selection Fee (ARUP test code 3002076) will be added to orders that utilize the first option. If multiple specimens are sent to ARUP without a request for pathologist block/slide selection or individual orders, they will be held until clarification is provided.

Stability: FFPE Tumor Tissue: Ambient: Indefinitely; Refrigerated: Indefinitely; Frozen: Unacceptable

Methodology: Quantitative Real-Time Polymerase Chain Reaction

Note:

CPT Codes: 81305

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Refer to report.

Reference Interval:

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.