

NEW TEST

[Click for Pricing](#)

Xylazine, Urine

3021300, XYLAZINE U

Specimen Requirements:

Patient Preparation:

Collect: Urine.

Specimen Preparation: Transfer 1 mL urine to an ARUP standard transport tube. (Min: 0.3 mL) Test is not performed at ARUP; separate specimens must be submitted when multiple tests are ordered.

Transport Temperature: Refrigerated. Also acceptable: Room temperature or frozen.

Unacceptable Conditions:

Remarks:

Stability: Ambient: 1 month; Refrigerated: 1 month; Frozen: 1 month

Methodology: Quantitative High Performance Liquid Chromatography-Tandem Mass Spectrometry

Note:

CPT Codes: 80375

New York DOH Approval Status: This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Refer to report

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.