

NEW TEST

[Click for Pricing](#)

Allergen, Food, Walnut (*Juglans* spp) Components IgE

3021309, WALNUT COM

Specimen Requirements:

Patient Preparation:	Multiple patient encounters should be avoided.
Collect:	Plain red or serum separator tube (SST)
Specimen Preparation:	Separate serum from cells ASAP or within 2 hours of collection. Transfer 1.0 mL serum to an ARUP standard transport tube. (Min: 0.75 mL). For multiple allergen orders refer to "Allergen Specimen Collection Instruction" at www.aruplab.com/testing/resources/specimen .
Transport Temperature:	Refrigerated.
Unacceptable Conditions:	Postmortem samples
Remarks:	
Stability:	After separation from cells: Ambient: 48 hours; Refrigerated: 2 weeks; Frozen: 1 month
Methodology:	Quantitative ImmunoCAP Fluorescent Enzyme Immunoassay
Note:	
CPT Codes:	86008x2
New York DOH Approval Status:	This test is New York DOH approved.

Interpretive Data:

Reference Interval:

Test Number	Components	Reference Interval
	Allergen, Walnut Component, rJug r 3	Less than or equal to 0.09 kU/L
	Allergen, Walnut rJug r 1	Less than or equal to 0.09 kU/L

HOTLINE NOTE: Refer to the Hotline Test Mix for interface build information.