



INFLUENZA, SARS-COV-2, AND RSV BY NAA SPECIMEN COLLECTION AND SHIPPING INSTRUCTIONS





Please note that information is changing rapidly, and ARUP will update this document frequently. For the most up-to-date version of this document, please visit <https://ltd.aruplab.com/Tests/Pub/3002638>

IMPORTANT NOTES:



- A completed [Patient Demographics Form for Public Health Reporting](#) is required. Failure to provide this information may result in a follow-up call from your state or local health department.
- **Specimen source is required.**

Test Overview: The ARUP Influenza, SARS-Cov-2, and RSV by NAA test is intended for the qualitative detection and differentiation of 2019 novel coronavirus SARS-CoV-2, influenza A/B virus, and/or RSV virus RNA from individuals suspected of respiratory viral infection consistent with COVID-19.

Reported in 1–4 days

Specimen Sources	Nasopharyngeal swab	Oropharyngeal swab	Nasal swab
VTM/UTM or equivalent* 	X	X	X
Saline 	X	X	X
Aptima Multitest Swab Specimen Collection Kit 		X	X
Aptima Specimen Transfer Kit 	X	X	X

INFLUENZA, SARS-COV-2, AND RSV BY NAA SPECIMEN COLLECTION AND SHIPPING INSTRUCTIONS

Specimen Sources	Nasopharyngeal swab	Oropharyngeal swab	Nasal swab
cobas PCR Media Uni & Dual Swab Sample Kits 			X
eSwab 	X	X	X

***VTM/UTM equivalent media: M4, M4-RT, M5, M5-RT, M6, UTM, BD Universal Transport Media (UVT), Puritan UniTranz media, Quest VCM (green lid), Starplex Multitrans System, UTM/VTM manufactured or aliquoted in the ARUP Reagent Laboratory, Hardy Diagnostics, VTM, Xpert Viral Transport Media**

Unacceptable media (including but not limited to): Media with guanidine-containing materials, 'molecular media' that indicates inactivation of pathogens and preservation of RNA/DNA, charcoal media: Amies gel swabs, Meridian Viral Transport Media, Beaver Nano-technologies cell-free DNA tube, PrimeStore MTM, CY-F003-40 viral transport media, DNA/RNA Shield, Abbott Molecular Multi-Collect Transport Tube, Ruhof 354CVTK, MANTACC UTM Transport Medium.

Unacceptable Specimens and Supplies: Swabs not in media; wood swabs; calcium alginate swabs; specimens sent in tubes with pop-top lids/caps; specimens in glass tubes.

If submitting media that is not listed in this guide, please send the package insert or ingredient list to COVID-19Requests@aruplab.com so it can be reviewed for compatibility with our testing platforms.

Specimen Type: Nasopharyngeal swab in a collection tube with viral transport medium (or equivalent medium). Also acceptable: oropharyngeal or nasal swab.

Nasopharyngeal swab instructions:

- Seat the patient comfortably and tilt the head back.
- Insert a nasal speculum, if available.
- Insert a swab through the nares until resistance is met due to contact with the nasopharynx.
- Rotate swab gently and maintain contact with the nasopharynx for 20–30 seconds or until coughing is induced.
- Break off the swab at the indicated spot, ensuring it fits completely within the collection tube.
- Close the tube, making sure the lid is secure.

Oropharyngeal swab instructions:

- Seat the patient comfortably and open the mouth.
- Insert a swab into the back of the throat.
- Rub the swab gently over the posterior oropharynx and both tonsillar pillars.
- Avoid touching other surfaces within the mouth, such as the tongue.
- Break off the swab at the indicated spot, ensuring it fits completely within the collection tube.
- Close the tube, making sure the lid is secure.

INFLUENZA, SARS-COV-2, AND RSV BY NAA SPECIMEN COLLECTION AND SHIPPING INSTRUCTIONS

Nasal swab instructions:

- Seat the patient comfortably and tilt the head back.
- Insert the swab into the first nostril until resistance is met at the level of the turbinates (approximately 1 inch into the nostril).
- Rotate the swab a few times against the nasal wall and remove from nostril.
- Using the same swab, carefully insert the swab into the second nostril until resistance is met at the level of the turbinates (approximately 1 inch into the nostril).
- Rotate the swab a few times against the nasal wall and remove from nostril.
- Break off the swab at the indicated spot, ensuring it fits completely within the collection tube.
- Close the tube, making sure the lid is secure.

IMPORTANT NOTES:

- Leaking tubes cannot be processed.
- Appropriate transport medium must be used.

Specimen Preparation: Please work with your laboratory management personnel, frontline clinicians, and collection locations to ensure the proper collection of Influenza, SARS-Cov-2, and RSV by NAA specimens. Specimens that have leaked in transport **cannot be processed** and will need to be recollected.

- Label the tube with at least two patient identifiers and place specimen in an individually sealed bag.
- Freeze and keep frozen.

When preparing Influenza, SARS-Cov-2, and RSV by NAA swab specimens, please make sure to snap off the shaft of the swab at the scored indicator line, leaving the flocked end of the swab inside the collection tube. Securely recap the collection tube. Please work with your collection staff to determine the appropriate procedure for cutting swabs that do not have scored indicator lines. The shaft of the swab should not be in contact with the lid, as this prevents a secure seal.

Do not place films, such as paraffin film (parafilm), between the tube and the lid or around the lid of the transport container. Failure to prepare specimens correctly may cause leakage and will prevent tubes from being processed.

Storage/Transport Temperature: Frozen

Specimen Stability: Swabs: Ambient: unacceptable; refrigerated: 2 days; frozen: 1 month

Shipping Instructions: For the most up-to-date transport instructions, refer to the [ARUP Laboratory Test Directory](#). The CDC states that specimens should be handled, packed, and shipped according to UN 3373 Biological Substance, Category B requirements, which is the standard process for shipping specimens to ARUP (refer to the [ARUP Specimen Transport Guide](#)). Current transportation practices have not changed. Additional information is available in the [CDC's biosafety guidelines](#).

How to Order: Refer to Influenza, SARS-Cov-2, and RSV by NAA (test code 3003443) in the test directory. The order must be placed by an authorized healthcare professional within the United States.

How to Obtain Supplies: When supply chains allow, supplies are available through eSupply using ARUP Connect™. Clients may obtain their own supplies. Acceptable equivalents for transport medium and swabs are listed above.

Information:

2019 Novel Coronavirus SARS-CoV-2 information from ARUP: aruplab.com/infectious-disease/coronavirus

For Client Assistance:

ARUP Client Services: 800-522-2787
