

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 7/31/1943
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Heavy Metals Panel 4, Urine with Reflex to Arsenic Fractionated

ARUP test code 0020572

Hours Collected	Random hr	
Total Volume	Random mL	
Creatinine, Urine - per volume	258 mg/dL	
Creatinine, Urine - per 24h	Not Applicable mg/d	(Ref Interval: 500-1400)
Cadmium, Urine - per volume	3.9 ug/L H	(Ref Interval: 0.0-1.0)
	INTERPRETATION INFORMATION: Cadmium, Urine Urine cadmium levels can be used to assess cadmium body burden. In chronic exposures, the kidneys are the primary target organ. Symptoms associated with cadmium toxicity vary based upon route of exposure and may include tubular proteinuria, fever, headache, dyspnea, chest pain, conjunctivitis, rhinitis, sore throat and cough. Ingestion of cadmium in high concentration may cause vomiting, diarrhea, salivation, cramps, and abdominal pain. See Compliance Statement B: aruplab.com/CS	
Cadmium, Urine - per 24h	Not Applicable ug/d	(Ref Interval: 0.0-3.2)
Cadmium, Urine - ratio to CRT	1.5 ug/g CRT	(Ref Interval: 0.0-3.2)
Lead, Urine - per volume	<5.0 ug/L	(Ref Interval: 0.0-5.0)
	INTERPRETIVE INFORMATION: Lead, Urine Quantification of urine excretion rates before or after chelation therapy has been used as an indicator of lead exposure. Urinary excretion of >125 mg of lead per 24 hours is usually associated with related evidence of lead toxicity. See Compliance Statement B: aruplab.com/CS	

H=High, L=Low, *=Abnormal, C=Critical

Lead, Urine - per 24h	Not Applicable ug/d	(Ref Interval: 0.0-8.1)
Lead, Urine - ratio to CRT	Not Applicable ug/g CRT	(Ref Interval: 0.0-5.0)
	Unable to accurately calculate the creatinine normalized result due to a low per volume result.	
Mercury, Urine - per volume	3.0 ug/L	(Ref Interval: 0.0-5.0)
	INTERPRETIVE INFORMATION: Mercury, Urine	
	Urinary mercury levels predominantly reflect acute or chronic elemental or inorganic mercury exposure. Urine concentrations in unexposed individuals are typically less than 10 ug/L. 24 hour urine concentrations of 30 to 100 ug/L may be associated with subclinical neuropsychiatric symptoms and tremors. Concentrations greater than 100 ug/L can be associated with overt neuropsychiatric disturbances and tremors. Urine mercury levels may be useful in monitoring chelation therapy.	
	Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS	
Mercury, Urine - per 24h	Not Applicable ug/d	(Ref Interval: 0.0-20.0)
Mercury, Urine - ratio to CRT	1.2 ug/g CRT	(Ref Interval: 0.0-20.0)
Arsenic Urine - per volume	37.0 ug/L H	(Ref Interval: 0.0-34.9)
	INTERPRETIVE INFORMATION: Arsenic, Urine w/ Reflex to Fractionated	
	The ACGIH Biological Exposure Index (BEI) for arsenic in urine is 35 ug/L. The ACGIH BEI is based on the sum of inorganic and methylated species. For specimens with a total arsenic concentration of 35 to 2000 ug/L, fractionation is automatically performed to determine the proportions of inorganic, methylated and organic species. It may be appropriate to request fractionation for specimens with total arsenic greater than 30 ug/gCRT despite a total arsenic concentration less than 35 ug/L. If low-level chronic poisoning is suspected, the ug/gCRT ratio may be a more sensitive indicator of arsenic exposure than the total arsenic concentration.	
	Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS	
Arsenic Urine - per 24h	Not Applicable ug/d	(Ref Interval: 0.0-49.9)
Arsenic, Urine - ratio to CRT	14.3 ug/g CRT	(Ref Interval: 0.0-29.9)

H=High, L=Low, *=Abnormal, C=Critical

Arsenic, Fractionated, Urine

ARUP test code 0020734

Arsenic, Organic	34.3 ug/L
------------------	-----------

Arsenic, Inorganic	<10.0 ug/L
--------------------	------------

Arsenic, Methylated	<10.0 ug/L
---------------------	------------

INTERPRETIVE INFORMATION: Arsenic, Fractionated Urine

The ACGIH Biological Exposure Index for the sum of inorganic and methylated species of arsenic is 35 ug/L. Inorganic species of arsenic are most toxic. Methylated species arise primarily from metabolism of inorganic species but may also come from dietary sources and are of moderate toxic potential. The organic species of arsenic are considered nontoxic and arise primarily from food. The sum of the inorganic, methylated, and organic species of arsenic may be lower than the total arsenic concentration due to the presence of unidentified organic species of arsenic.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS

H=High, L=Low, *=Abnormal, C=Critical

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Hours Collected	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Total Volume	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Creatinine, Urine - per volume	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Creatinine, Urine - per 24h	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Cadmium, Urine - per volume	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Cadmium, Urine - per 24h	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Cadmium, Urine - ratio to CRT	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Lead, Urine - per volume	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Lead, Urine - per 24h	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Lead, Urine - ratio to CRT	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Mercury, Urine - per volume	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Mercury, Urine - per 24h	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Mercury, Urine - ratio to CRT	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Arsenic Urine - per volume	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Arsenic Urine - per 24h	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Arsenic, Urine - ratio to CRT	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Arsenic, Organic	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Arsenic, Inorganic	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Arsenic, Methylated	19-142-127239	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

ARUP LABORATORIES | 800-522-2787 | aruplab.com
500 Chipeta Way, Salt Lake City, UT 84108-1221
Tracy I. George, MD, Laboratory Director

Patient: Patient, Example
ARUP Accession: 19-142-127239
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Page 4 of 4 | Printed: 2/19/2021 2:57:45 PM
4848