

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 10/10/1979
Gender: Male
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Nickel, Urine

ARUP test code 0025045

| | | |
|--------------------------------|---|---------------------------|
| Hours Collected | 24 hr | |
| | Per 24h calculations are provided to aid interpretation for collections with a duration of 24 hours and an average daily urine volume. For specimens with notable deviations in collection time or volume, ratios of analytes to a corresponding urine creatinine concentration may assist in result interpretation. | |
| Total Volume | 1500 mL | |
| Creatinine, Urine - per volume | 100 mg/dL | |
| Creatinine, Urine - per 24h | 1500 mg/d | (Ref Interval: 1000-2500) |
| Nickel, Urine - per volume | 5.0 ug/L | (Ref Interval: 0.0-10.4) |
| | TEST INFORMATION: Nickel, Urine Measurement of nickel is not recommended in asymptomatic individuals or in individuals with a low likelihood of exposure. Elevations in nickel urine should be interpreted with caution in individuals with no exposure risks, and may indicate contamination of the specimen. This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes. | |
| Nickel, Urine - per 24h | 7.5 ug/d | (Ref Interval: 0.0-14.9) |
| Nickel, Urine - ratio to CRT | 5.0 ug/g CRT | (Ref Interval: 0.0-9.9) |

H=High, L=Low, *=Abnormal, C=Critical

| VERIFIED/REPORTED DATES | | | | |
|--------------------------------|---------------|------------------|------------------|-------------------|
| Procedure | Accession | Collected | Received | Verified/Reported |
| Hours Collected | 22-221-112499 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00 |
| Total Volume | 22-221-112499 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00 |
| Creatinine, Urine - per volume | 22-221-112499 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00 |
| Creatinine, Urine - per 24h | 22-221-112499 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00 |
| Nickel, Urine - per volume | 22-221-112499 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00 |
| Nickel, Urine - per 24h | 22-221-112499 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00 |
| Nickel, Urine - ratio to CRT | 22-221-112499 | 00/00/0000 00:00 | 00/00/0000 00:00 | 00/00/0000 00:00 |

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at:

ARUP LABORATORIES | 800-522-2787 | aruplab.com
500 Chipeta Way, Salt Lake City, UT 84108-1221
Jonathan R. Genzen, MD, PhD, Laboratory Director

Patient: Patient, Example
ARUP Accession: 22-221-112499
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Page 2 of 2 | Printed: 8/16/2022 8:48:40 AM
4848