

Client: Example Client ABC123  
123 Test Drive  
Salt Lake City, UT 84108  
UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example**

**DOB:** 6/16/1991  
**Gender:** Male  
**Patient Identifiers:** 01234567890ABCD, 012345  
**Visit Number (FIN):** 01234567890ABCD  
**Collection Date:** 00/00/0000 00:00

**B-Type Natriuretic Peptide**

ARUP test code 0030191

B-Type Natriuretic Peptide

15 pg/mL (Ref Interval: <=100)

**INTERPRETIVE INFORMATION:** B-Type Natriuretic Peptide  
A cutoff of 100 pg/mL has been demonstrated to provide the maximal combination of sensitivity, specificity and negative predictive value for contributing to the diagnosis of congestive heart failure (CHF). A BNP value greater than 100 pg/mL is consistent with a diagnosis of CHF in the appropriate clinical setting. False positive results are more common in females greater than 75 years of age. Blood concentrations of natriuretic peptides may also be elevated in patients with myocardial infarction and in patients that are candidates for or are undergoing renal dialysis.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
B-Type Natriuretic Peptide	23-109-104834	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: