

Client: Example Client ABC123 123 Test Drive Salt Lake City, UT 84108 UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB 1/4/1961 **Gender:** Female

Patient Identifiers: 01234567890ABCD, 012345

Visit Number (FIN): 01234567890ABCD **Collection Date:** 00/00/0000 00:00

Striated Muscle Antibodies, IgG with Reflex to Titer

ARUP test code 0050746

Striated Muscle Antibodies, IgG Screen

Detected * (Ref Interval: <1:40)

Striated Muscle Antibodies, IgG detected. Titer results to follow.

INTERPRETIVE DATA: Striated Muscle Antibodies, IgG Screen

In the presence of acetylcholine receptor (AChR) antibody, striated muscle antibodies, which bind in a cross-striational pattern to skeletal and heart muscle tissue sections, are associated with late-onset myasthenia gravis (MG). Striated muscle antibodies recognize epitopes on three major muscle proteins, including: titin, ryanodine receptor (RyR) and Kv1.4 (an alpha subunit of voltage-gated potassium channel [VGKC]). Isolated cases of striated muscle antibodies may be seen in patients with certain autoimmune diseases, rheumatic fever, myocardial infarction, and following some cardiotomy procedures.

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

Striated Muscle Antibodies, IgG Titer (Reflex for 0050746 STM R - Do NOT give this test code to clients - Panel Component ONLY.)

ARUP test code 2012516

Striated Muscle Antibodies, IgG Titer

1:320 * (Ref Interval: <1:40)

INTREPRETIVE INFORMATION: Striated Muscle Abs, IgG Titer

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

H=High, L=Low, *=Abnormal, C=Critical

4848



VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
Striated Muscle Antibodies, IgG Screen	22-189-144237	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Striated Muscle Antibodies, IgG Titer	22-189-144237	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical