

Client: Example Client ABC123 123 Test Drive Salt Lake City, UT 84108 UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example** 

DOB 2/28/1963
Gender: Female

**Patient Identifiers:** 01234567890ABCD, 012345

**Visit Number (FIN):** 01234567890ABCD **Collection Date:** 00/00/0000 00:00

## Rheumatoid Factors, IgA, IgG, and IgM by ELISA

ARUP test code 0051298

Rheumatoid Factor, IgA by ELISA 12 Units H (Ref Interval: <=6)

INTERPRETIVE INFORMATION: Rheumatoid Factor, IgA by ELISA

The presence of all three rheumatoid factor (RF) isotypes at abnormal levels has high specificity for a diagnosis of rheumatoid arthritis (RA). However, the presence of RF isotypes in any combination may be found in a variety of conditions, including Sjogren syndrome and hepatitis infections.

Rheumatoid Factor, IgG by ELISA 7 Units H (Ref Interval: <=6)

INTERPRETIVE INFORMATION: Rheumatoid Factor, IgG by ELISA

The presence of all three rheumatoid factor (RF) isotypes at abnormal levels has high specificity for a diagnosis of rheumatoid arthritis (RA). However, the presence of RF isotypes in any combination may be found in a variety of conditions, including Sjogren syndrome and hepatitis infections.

Rheumatoid Factor, IgM by ELISA 12 Units H (Ref Interval: <=6)

INTERPRETIVE INFORMATION: Rheumatoid Factor, Igm by ELISA

The presence of all three rheumatoid factor (RF) isotypes at abnormal levels has high specificity for a diagnosis of rheumatoid arthritis (RA). However, the presence of RF isotypes in any combination may be found in a variety of conditions, including Sjogren syndrome and hepatitis infections.

H=High, L=Low, \*=Abnormal, C=Critical

4848



VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
Rheumatoid Factor, IgA by ELISA	22-189-400909	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Rheumatoid Factor, IgG by ELISA	22-189-400909	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Rheumatoid Factor, IgM by ELISA	22-189-400909	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical