

Client: Example Client ABC123
 123 Test Drive
 Salt Lake City, UT 84108
 UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 1/21/1959
Gender: Male
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Cyclosporine A, 2-Hour Post Dose (C2) by Tandem Mass Spectrometry

ARUP test code 0058902

CyclosporineA, 2-Hr Post(C2), HPLC-MS/MS 742.0 ng/mL

INTERPRETIVE INFORMATION: Cyclosporine A, 2-Hr Post(C2)
 HPLC-MS/MS

Cyclosporine A levels in specimens drawn 2 hours post-dose may estimate the AUC better than trough specimens. The optimal therapeutic range for a given patient may differ from this suggested range based on the indication for therapy, treatment phase (initiation or maintenance), use in combination with other drugs, time of specimen collection relative to prior dose, type of transplanted organ, and/or the therapeutic approach of the transplant center. A suggested target range for renal transplant is 800-1700 ng/mL. A suggested target range for liver transplant is 600-1000 ng/mL. The higher ranges represent concentrations immediately post-transplant and the lower ranges represent concentrations during the maintenance phase.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
CyclosporineA, 2-Hr Post(C2), HPLC-MS/MS	19-273-115035	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical