

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 7/8/1973
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Insulin, Other

ARUP test code 0070022

Insulin, Other

5 uIU/mL

INTERPRETIVE INFORMATION: Insulin, Other

This test reacts on a nearly equimolar basis with the analogs insulin aspart, insulin glargine, and insulin lispro. Insulin detemir exhibits approximately 50 percent cross-reactivity. Test reactivity with insulin glulisine is negligible (< 3 percent). To convert to pmol/L, multiply uIU/mL by 6.0. The reference interval for fasting insulin is 3-25 uIU/mL

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Insulin, Other	23-305-126934	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical