

Client: Example Client ABC123 123 Test Drive

Salt Lake City, UT 84108 UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example** 

**DOB** 7/2/1979 **Gender:** Female

**Patient Identifiers:** 01234567890ABCD, 012345

**Visit Number (FIN):** 01234567890ABCD **Collection Date:** 00/00/0000 00:00

## **Insulin, 30 Minutes**

ARUP test code 0070064

Insulin, 30 Minutes

120 uIU/mL H (Ref Interval: 20-112)

INTERPRETIVE INFORMATION: Insulin 30 Minutes This test reacts on a nearly equimolar basis with the analogs insulin aspart, insulin glargine, and insulin lispro. Insulin detemir exhibits approximately 50 percent cross-reactivity. Test reactivity with insulin glulisine is negligible (<3 percent). To convert to pmol/L, multiply uIU/mL by 6.0. The reference interval is based on a 75 g. glucose challenge.

VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
Insulin, 30 Minutes	23-311-145799	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical

4848