

Client: Example Client ABC123  
123 Test Drive  
Salt Lake City, UT 84108  
UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example**

**DOB:** 7/24/2010  
**Gender:** Female  
**Patient Identifiers:** 01234567890ABCD, 012345  
**Visit Number (FIN):** 01234567890ABCD  
**Collection Date:** 01/01/2017 12:34

**Insulin, 30 Minutes**

ARUP test code 0070064

Insulin, 30 Minutes

**467 uIU/mL H (Ref Interval: 20-112)**

INTERPRETIVE INFORMATION: Insulin 30 Minutes  
This test reacts on a nearly equimolar basis with the analogs insulin aspart, insulin glargine, and insulin lispro. Insulin detemir exhibits approximately 50 percent cross-reactivity. Test reactivity with insulin glulisine is negligible (<3 percent). To convert to pmol/L, multiply uIU/mL by 6.0. The reference interval is based on a 75 g. glucose challenge.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Insulin, 30 Minutes	19-309-102235	11/5/2019 9:33:00 AM	11/6/2019 5:53:53 PM	11/7/2019 12:12:00 AM

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: