

Client: Example Client ABC123

123 Test Drive

Salt Lake City, UT 84108 UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB 2/1/1981
Gender: Female

Patient Identifiers: 01234567890ABCD, 012345

Visit Number (FIN): 01234567890ABCD **Collection Date:** 00/00/0000 00:00

Insulin, 30 Minutes

ARUP test code 0070064

Insulin, 30 Minutes

104 uIU/mL

(Ref Interval: 20-112)

INTERPRETIVE INFORMATION: Insulin 30 Minutes
This test reacts on a nearly equimolar basis with the analogs
insulin aspart, insulin glargine, and insulin lispro. Insulin
detemir exhibits approximately 50 percent cross-reactivity.
Test reactivity with insulin glulisine is negligible (<3
percent). To convert to pmol/L, multiply uIU/mL by 6.0. The
reference interval is based on a 75 g. glucose challenge.

VERIFIED/REPORTED DATES				
Procedure	Accession	Collected	Received	Verified/Reported
Insulin, 30 Minutes	23-311-145815	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

4848