

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 6/3/2003
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Insulin, Free and Total

ARUP test code 0070155

Insulin Free **28 uIU/mL H** (Ref Interval: 3-25)

Total Insulin **36 uIU/mL H** (Ref Interval: 3-25)

INTERPRETIVE INFORMATION: Insulin, Free and Total
This test reacts on a nearly equimolar basis with the analogs insulin aspart, insulin glargine, and insulin lispro. Insulin detemir exhibits approximately 50 percent cross-reactivity. Test reactivity with insulin glulisine is negligible (<3 percent). To convert to pmol/L, multiply uIU/mL by 6.0. Reference intervals established for fasting specimens.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Insulin Free	21-305-145494	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Total Insulin	21-305-145494	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: