

Client: Example Client ABC123  
123 Test Drive  
Salt Lake City, UT 84108  
UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example**

**DOB:** 3/15/1971  
**Gender:** Female  
**Patient Identifiers:** 01234567890ABCD, 012345  
**Visit Number (FIN):** 01234567890ABCD  
**Collection Date:** 00/00/0000 00:00

**Insulin, Free and Total**

ARUP test code 0070155

Insulin Free 6 uIU/mL (Ref Interval: 3-25)

Total Insulin 7 uIU/mL (Ref Interval: 3-25)

**INTERPRETIVE INFORMATION: Insulin, Free and Total**  
This test reacts on a nearly equimolar basis with the analogs insulin aspart, insulin glargine, and insulin lispro. Insulin detemir exhibits approximately 50 percent cross-reactivity. Test reactivity with insulin glulisine is negligible (<3 percent). To convert to pmol/L, multiply uIU/mL by 6.0. Reference intervals established for fasting specimens.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Insulin Free	23-315-108743	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Total Insulin	23-315-108743	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: