

Client: Example Client ABC123
123 Test Drive
Salt Lake City, UT 84108
UNITED STATES

Physician: Doctor, Example

Patient: Patient, Example

DOB: 8/8/1946
Gender: Female
Patient Identifiers: 01234567890ABCD, 012345
Visit Number (FIN): 01234567890ABCD
Collection Date: 00/00/0000 00:00

Thyroxine, Free by Equilibrium Dialysis/HPLC-Tandem Mass Spectrometry

ARUP test code 0093244

Free T4 by Equil Dialysis-TMS **3.7 ng/dL H** (Ref Interval: 1.1-2.4)

FREE T4 BY EQUIL DIALYSIS-TMS:
PREGNANCY REFERENCE INTERVALS

1ST TRIMESTER 0.7 - 2.0 ng/dL
2ND TRIMESTER 0.7 - 2.1 ng/dL
3RD TRIMESTER 0.5 - 1.6 ng/dL

INTERPRETIVE INFORMATION: FT4 ED-TMS

Some medications may induce transient changes in FT4 concentrations. This test is not recommended for patients currently on heparin treatment as FT4 concentrations may be falsely elevated.

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Free T4 by Equil Dialysis-TMS	24-004-100587	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, *=Abnormal, C=Critical

Unless otherwise indicated, testing performed at: