

Client: Example Client ABC123  
123 Test Drive  
Salt Lake City, UT 84108  
UNITED STATES

Physician: Doctor, Example

**Patient: Patient, Example**

**DOB:** 9/10/1941  
**Gender:** Female  
**Patient Identifiers:** 01234567890ABCD, 012345  
**Visit Number (FIN):** 01234567890ABCD  
**Collection Date:** 00/00/0000 00:00

**Immunoglobulins, Serum Cryoprecipitins (Reflex for 2002063 CRYGB QNT and 2002403 CRYO TYPING only. NOT orderable by clients.)**

ARUP test code 2002125

Immunoglobulin G, Cryoprecipitate **7 mg/dL H** (Ref Interval: 0-0)  
Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS

Immunoglobulin A, Cryoprecipitate **1 mg/dL H** (Ref Interval: 0-0)

Immunoglobulin M, Cryoprecipitate **26 mg/dL H** (Ref Interval: 0-0)

**Immunofixation Electrophoresis, Qualitative, Gel**

ARUP test code 2012572

Immunofix Electrophoresis Gel See Note  
The cryoglobulin is characterized as a Type II cryoglobulin with monoclonal IgM lambda and polyclonal IgG with an additional faint band in IgM kappa.  
INTERPRETIVE INFORMATION: Immunofixation Electrophoresis Gel  
This information should be correlated with the results of serum protein electrophoresis, quantitative immunoglobulins and other clinical and laboratory information.

EER Immunofix Electrophoresis Gel See Note  
Access ARUP Enhanced Report using either link below:  
-Direct access: [REDACTED]  
-Enter Username, Password: [REDACTED]  
Username: [REDACTED]  
Password: [REDACTED]

H=High, L=Low, \*=Abnormal, C=Critical

**Cryoglobulin, Qualitative with Reflex to IFE Typing and Quantitative IgA, IgG, and IgM**

ARUP test code 2002403

Cryoglobulin, Qualitative

**POS 24Hour \* (Ref Interval: NEG 72Hour)**  
 Cryoglobulin detected at 24 hrs.

Test developed and characteristics determined by ARUP Laboratories. See Compliance Statement B: aruplab.com/CS

VERIFIED/REPORTED DATES

Procedure	Accession	Collected	Received	Verified/Reported
Immunoglobulin G, Cryoprecipitate	19-178-101506	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Immunoglobulin A, Cryoprecipitate	19-178-101506	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Immunoglobulin M, Cryoprecipitate	19-178-101506	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Immunofix Electrophoresis Gel	19-178-101506	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
Cryoglobulin, Qualitative	19-178-101506	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00
EER Immunofix Electrophoresis Gel	19-178-101506	00/00/0000 00:00	00/00/0000 00:00	00/00/0000 00:00

END OF CHART

H=High, L=Low, \*=Abnormal, C=Critical